

Disinfectant

Helimatic® Disinfectant ... for automatic reprocessing of flexible endoscopes & medical devices

Properties

- For automated disinfection of flexible endoscopes and heat sensitive medical devices
- Broad efficacy spectrum
- Effective and economic
- Has a high level of material compatibility especially with plastics and metallic surfaces

Microbiological efficacy

Microbiocidal efficacy under clean conditions

Test organism / Virus	Conc.	Time	Temp.
Bactericidal (EN 13727, EN 14561)	1.0 %	5 min.	55 °C
Fungicidal (EN 13624, EN 14562)	1.0 %	5 min.	55 °C
Mycobactericidal (EN 14348, EN 14563)	1.0 %	5 min.	55 °C
Virucidal (EN 14476, 17111)	1.0 %	5 min.	55 °C

Expert reports are available on demand.

Reprocessing program proposal for flexible endoscopes: see page 25.

Product size	Ref.
5 l canister	18562
10 l canister	180349
200 l barrel	180350

Physico-Chemical Data

pH-value (20 °C):
Density (20 °C, g/cm³):
Appearance

Concentrate

4
Approx. 1.07
clear, liquid

Ready-to-use solution

4
Approx. 1.0

Dosage and instructions for use

The cleaning step plays a key role in thermo chemical processing. The success of the processing depends on thorough cleaning in order that the flexible endoscope is free of secretion- and protein-residuals. This prevents the incrustation of any remaining contamination during the subsequent disinfecting stage and drying phase.

Reprocessing of flexible endoscopes

- Wear gloves and personal protective equipment, follow the reprocessing recommendations of the endoscope manufacturer
- Pre cleaning in the examination room: immediately after the examination (enzymatic cleaner e.g. Helizyme)
- Manual cleaning in the reprocessing room: clean the channels and other parts of the endoscope with special cleaning brushes (Helizyme, enzymatic cleaner).
- Rinsing: Rinse with water
- Machine reprocessing in a washer-disinfector for flexible endoscopes.
- Allow to dry completely.
- Low temperature sterilization: if available and required
- Store in a cabinet recommended by the endoscope manufacturer.



Highly effective

Helimatic® Disinfectant– Composition:

100 g solution contains 20 g glutaral, solvents, salts of organic acids, corrosion inhibitors, excipients. Ingredients in accordance with the Regulations for Detergents EG 648/2004, < 5% phosphonate | Labelling of dangerous goods: see material safety data sheet (MSDS). Cautions: Use disinfectants safely. Always read the label and the product information before use. Do not use the product after the expiry date.