

Alkaline Detergent



Product Features:

- Removes soils and bioburden from stainless steel instruments, rubber, and plastics
- Accommodates instrument IFUs requiring the use of alkaline pH cleaners

Key Formula Ingredients:

- Potassium Hydroxide
- Potassium Silicate

Key Chemical Attributes

pH Neat	12.5-13.4
pH 1% Dilution	9.1-9.9
Density	Approximately 1.027
Color	Light Yellow

Handling and Storage:

Use with adequate ventilation. Avoid breathing vapor spray or mist. Store in a cool dry place. Keep partially used containers closed and away from heat and humidity. Check periodically for corrosion and integrity. Incompatible with acids, bases and oxidizers. Do not expose to extreme hot or cold temperatures. Keep out of reach of children. Keep in properly labeled containers.

Alkaline Cleaning Solution for Automated Washer/Disinfectors – Concentrate

Formulated for stainless steel surfaces, rubber, and plastic based items.

Safety Information:

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.




DANGER: May be corrosive to metals. Causes severe skin burns and eye damage. May cause damage to organs through prolonged or repeated exposure. Do not breath mist, vapor or spray. Wash hands and exposed areas thoroughly after handling. Wear eye and face protection, protective gloves and protective clothing.

Instructions for Use:

1. Automatically feed 1/8 to 2 oz of detergent per gallon of water (1 to 16 ml per liter of water)
2. For best results, water temperature should be 86°C to 130°F (30°C to 60°C)
3. Do not mix with other detergents
4. Follow facility policies and procedures

Product Ordering Information:

Product Description/ Name	Color Code	Order #	Item #	Pkg	Units Per Case
Alkaline Detergent		110-1530	2191002	2.5 gal	2
		110-1531	2191005	5 gal	1



For more information, visit belimed.com

REV7.5.2022_Q0

Engineers of Confidence.

Belimed
Infection Control

Belimed, Inc.
1.800.451.4118 | 1.843.216.7424
belimed.com