

Enzymatic Pretreatment Foam Spray



Specially formulated multi-enzyme dual foam/liquid pretreatment spray

For pre-cleaning and transportation of surgical grade instruments, reusable items, scopes and other devices

Product Features:

- · pH Neutral
- · Gentle on instruments
- · Easily rinsed
- · Fragrance-free
- · Stays wet up to 72 hours
- · Safe for end users
- · Reduces risk of biofilm formation

Key Formula Ingredients:

- Three (3) enzymes (protease for proteins, lipase for fatty acids and amylase for starches)
- $\cdot\;$ Unique blend of surfactants solubilize bioburden and retain wetness
- · Corrosion inhibitor

Key Chemical Attributes	
pH Neat	6.5-8.5
pH 1% Dilution	N/A
Density	1.010-1.020
Color	Blue

Material Compatibility:

The Belimed Protect® Enzymatic Pretreatment Foam Spray demonstrates wide range of material compatibility and is suitable for the following materials:

- Stainless Steel
- · Brass
- · Aluminum
- All commonly used plastics and elastomers, including Polyurethane, Silicone, and Teflon
- Anodized AluminumTitanium
- · Latexes and Synthetics (Vinyl and Nitrile)
- · Copper

Handling and Storage:

Do not store at extreme hot or cold temperatures.

Safety Information:

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.

Recommended Practice:

- · After each case, staff should spray all instrument trays
- · The instruments should be covered completely with product
- · After a few hours, check trays and re-spray if necessary

Instructions for Use:

- 1. Shake bottle before each use
- 2. Spray instruments and equipment as soon as possible after procedure
- Keep instruments and equipment wetted with solution until decontamination process is started
- 4. Follow facility policies and procedures

Product Ordering Information:

Product Description/ Name	Color Code	Order #	Item #	Pkg	Units Per Case
Enzymatic Pretreatment Foam Spray		106-8088	2110032	32 oz	12
		1008877	2111002	2 oz	48
		110-2632	2110001	1 gal	4
		110-2633	2110002	2.5 gal	2



For more information, visit belimed.com

REV10.8.2025_Q0

Engineers of Confidence.

