

Instrument Lubricant



Ready to use instrument lubricant

For lubrication and preparation of surgical instrument hinges and joints.

Product Features:

- Special formulation designed to lubricate and protect instruments, hinges, box locks, and general equipment
- Gentle, non-corrosive ingredients will not damage instruments
- Fragrance free
- Biodegradable
- pH neutral

Key Formula Ingredients:

- Food-grade mineral oil
- Unique blend of surfactants

Key Chemical Attributes	
pH Neat	5.5-7.5
pH 1% Dilution	6.0-7.0
Density	0.985-1.005
Color	Milky White

Material Compatibility:

The Belimed Protect® Instrument Lubricant demonstrates a wide range of material compatibility and is suitable for the following materials:

- Stainless Steel
- Aluminum
- Anodized Aluminum
- Titanium
- Copper
- Brass
- All commonly used plastics and elastomers, including Polyurethane, Silicone, and Teflon
- Latexes and Synthetics (Vinyl and Nitrile)

Handling and Storage:

Do not store at extreme hot or cold temperatures.


Safety Information:

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.

Instructions for Use:

1. Shake bottle before each use
2. Apply evenly by holding spray nozzle 6-8 inches (15-20 cm) over surgical instruments. Rinsing is not required
3. Follow facility policies and procedures and instrument Instructions for Use

Product Ordering Information:

Product Description/ Name	Color Code	Order #	Item #	Pkg	Units Per Case
Instrument Lubricant		106-8091	2152032	32 oz	12



For more information, visit belimed.com

REV12.22.2020_Q0

Engineers of Confidence.