

Ophthalmic Neutral Detergent



The first and only neutral pH non-enzymatic detergent formulated specifically for ophthalmic instrumentation

For manual or automated cleaning of surgical instruments and other items used for eye surgery.

Product Features:

- Made specifically for ophthalmic instrumentation based on OICS Task Force guidelines
- Independently tested and classified as non-irritating to skin and eyes and non-cytotoxic in its undiluted form
- pH neutral
- Non-corrosive
- Gentle on instruments
- Fragrance free
- Biodegradable and aquatic life safe

Key Formula Ingredients:

- Unique blend of surfactants to solubilize bioburden
- Corrosion inhibitor
- Sequestering agent to combat hard water minerals

Key Chemical Attributes	
pH Neat	6.5-8.0
pH 1% Dilution	7.0-7.5
Density	1.020-1.030
Color	Straw

Material Compatibility:

The Belimed Protect[®] Ophthalmic Neutral Detergent demonstrates a wide range of material compatibility and is suitable for the following materials:

- Stainless Steel
- Aluminum
- Anodized Aluminum
- Titanium
- Copper
- Brass
- All commonly used plastics and elastomers, including Polyurethane, Silicone, and Teflon
- Latexes and Synthetics (Vinyl and Nitrile)

Handling and Storage:

Do not store at extreme hot or cold temperatures.




Safety Information:

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.

Instructions for Use - Manual/Ultrasonic/Washer Disinfector:

1. Add 1/8 to 1 oz of detergent per gallon (1 to 8 ml/L) of water
2. For best results, ensure water temperature is 80°-140°F (26°-60°C)
3. Process instruments according to facility policies and procedures and the instrument manufacturer IFU

Product Ordering Information:

Product Description/Name	Color Code	Order #	Item #	Pkg	Units Per Case
Ophthalmic Neutral Detergent		102-0055	2193001	1 gal	4
		102-0056	2193002	2 gal	2
		102-0057	2193005	5 gal	1



For more information, visit belimed.com

REV12.22.2020_Q0

Engineers of Confidence.