



## 565 PLUS Ultrasonic Series

The 565 PLUS Ultrasonics are easy to use systems, for cleaning and rinsing a variety of instruments - from heavy orthopedic sets to cannulated instruments and daVinci endowrists. The PLUS has a two tray capacity, while the PLUS3 has a three tray capacity.

Engineers of Confidence.



# 565 PLUS Ultrasonic Series

### Accurate Detergent Dosing

- Accurate, repeatable flow meter dosing
- Detergent filter before flowmeter
- Visual detergent level indicator on display
- Reset button for detergent level

### Effective Cleaning

- Pulsed flushing pump for improved cleaning
- Improved cannulated flushing algorithm
- Integrated water mixing valve
- Set water temperature from display
- Additional ultrasonic transducers and power with the PLUS3, 3 tray model

### Reliable and User Friendly


- Electronic water level sensors
- Intuitive display flow and information
- Data storage USB port
- Password protected access levels
- Optional pump for wall mounted drain
- Powered lid with 70 lb tray lift capacity

### Safety Glass Lid

- Visually confirm instruments are submerged
- Allows manual lifting of the lid
- Damped closing feature when released

### Instrument Connection Options

- Press fit for daVinci endowrists and cannulated devices
- Twist luer lock
- Flexible diaphragm for non-luer, cannulated devices



The PLUS has a two tray capacity, while the PLUS3 has a three tray capacity.



### Standard Tray

- Two trays per load for PLUS
- Three trays per load for PLUS3



### Cannulated Instrument Tray

- 10 connections per tray
- Wash and rinse flushing
- Two trays per load for PLUS
- Two trays plus additional standard tray per load for PLUS3



### daVinci Tray

- 10 instrument capacity per tray
- Wash and rinse flushing
- Two trays per load for PLUS
- Two trays plus additional standard tray per load for PLUS3

REV9.28.2020\_Q320

Engineers of Confidence.

**Belimed**  
Infection Control

Belimed, Inc.  
1.800.451.4118 | 1.843.216.7424  
[belimed.com](http://belimed.com)