

Cart Rinse Plus



Cart Rinse Plus for Cart Washers and Automated Washer/Disinfectors

Concentrated cart rinse aid for spot free case carts and stainless steel products.

Product Features:

- Ensures complete drying of case carts in automated cart washers and other automated equipment
- Conditions water for spot free drying and sheeting action reduces dry time
- Concentrated formula optimizes dosing levels
- Gentle ingredients will not damage case carts, casters or containers
- Biodegradable and aquatic life safe
- Fragrance free

Handling and Storage:

Do not store at extreme hot or cold temperatures.

Safety Information:

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.

Instructions for Use:

1. Add 1/16 to 2/5 oz of rinse aid (0.5 to 3 ml/L) per gallon of water.
2. Follow facility policies and procedures.

Key Chemical Attributes

pH Neat	2.0-3.5
pH Use Dilution	5-7
Density	Approximately 1.0075
Color	Blue

Product Ordering Information:

Product Description/ Name	Color Code	Order #	Item #	Pkg	Units Per Case
Cart Rinse Plus		111-1905	2171002	2.5 gal	2
		111-1909	2171005	5 gal	1
		111-1910	2171015	15 gal	1



For more information, visit belimed.com

REV7.5.2022_Q0

Engineers of Confidence.