

Neutralizing Detergent



Neutralizing Detergent for Sink, Ultrasonic and Automated Washer Disinfector – Concentrate

For neutralizing alkaline detergent residues and improving rinse water quality.

Product Features:

- Neutralizes alkaline detergent residues
- Removes dissolved mineral content from tap water rinse
- Eliminates white, powdery deposits from instrument surfaces, rubber and plastic items
- Rejuvenates surfaces of surgical instruments

Key Formula Ingredients:

- Citric acid

Key Chemical Attributes

pH Neat	2.1-3.1
pH 1% Dilution	N/A
Density	1.05-1.06
Color	Clear

Handling and Storage:

Use with adequate ventilation. Avoid breathing vapor spray or mist. Keep partially used containers closed and away from heat and humidity. Do not store in extreme hot or cold temperatures. Keep out of reach of children.

Safety Information:

Before using product please refer to Instructions for Use, Safety Data Sheet and product label.




DANGER: Causes serious eye irritation. Wash skin thoroughly after handling. Wear eye/face protection, protective gloves and protective clothing.

Handling and Storage:

1. Automatically feed 1/8 to 2 oz of detergent per gallon of water (1 to 16 ml per liter of water).
2. For best results, water temperature should be 86°F to 140°F (30°C to 60°C).
3. Follow facility policies and procedures.

Product Ordering Information:

Product Description/ Name	Color Code	Order #	Item #	Pkg	Units Per Case
Neutralizing Detergent		110-1533	2192002	2.5 gal	2
		110-1534	2192005	5 gal	1



For more information, visit belimed.com

REV12.22.2020_Q0

Engineers of Confidence.