

### Neutral Non-enzymatic Pretreatment Spray

#### SECTION 1: PRODUCT AND COMPANY IDENTIFICATION

Product identifier used on the label : **Neutral Non-enzymatic Pretreatment Spray**

Other means of identification : Instrument pretreatment foam/liquid spray for surgical grade instruments, scopes and medical devices  
: 2110016, 1032459

#### Recommended use of the chemical and restrictions on use:

Product for hospital and professional use only for the cleaning of surgical instruments. Not for home use.

#### Name, address, and telephone number of the chemical manufacturer:

BELIMED INC.  
8351 Palmetto Commerce Parkway, Suite 101  
Ladson, SC 29456 USA  
Telephone Number for Information : 1-800-451-4118

#### Emergency telephone number

CHEMTREC : 1-800-424-9300 (24/7, US and Canada)

#### SECTION 2: HAZARD IDENTIFICATION

**Warning** : Product is not considered hazardous or irritating.

**Classification of the substance or mixture** : Not a hazardous substance or mixture.

#### Signal word, hazard statement(s), symbol(s) and precautionary statement(s)

**Signal word** : None  
Product is not considered hazardous or irritating according to OSHA regulations. However, personal protective equipment should be worn by OSHA bloodborne pathogen guidelines when handling soiled instruments

**Precautionary Statements** : Prevention: Wash hands thoroughly after handling.  
: Response: Get medical advice/ attention if you feel unwell.

**Storage** : Store in accordance with local regulations.

**Disposal** : Dispose of contents/container to an approved waste disposal facility.

**Hazards not otherwise classified** : Not known.

**Percentage of ingredient(s) with unknown toxicity** : None known.

Classification complies with OSHA Hazard Communication Standard (29 CFR 1910.1200) and is consistent with the provisions of the United Nations Globally Harmonized System of Classification and Labeling of Chemicals (GHS).

## Neutral Non-enzymatic Pretreatment Spray

### SECTION 3: COMPOSITION / INFORMATION ON INGREDIENTS

Pure substance/mixture : Mixture

### SECTION 4: FIRST AID MEASURES

#### Description of necessary measures

**Inhalation:** Remove from exposure area to fresh air. Seek medical attention immediately.

**Skin contact:** Immediately wash skin with soap and water for at least 20 minutes while removing contaminated clothing. Seek medical attention if irritation persists.

**Eye contact:** Immediately flush with cool water for at least 20 minutes. Seek medical attention if irritation persists.

**Ingestion:** Do NOT induce vomiting. Dilute by rinsing the mouth and giving cool water to drink. Never give anything by mouth to an unconscious person. Contact physician or local poison control center.

#### Most important symptoms and effects, both acute and delayed

**Inhalation:** May cause allergic respiratory reaction.

**Skin contact:** Health injuries are not known or expected.

**Eye contact:** Health injuries are not known or expected.

**Ingestion:** Ingestion may cause irritation of mouth, throat, and digestive tract, diarrhea, and vomiting.

#### Indication of any immediate medical attention and special treatment needed

**After inhalation:** Remove from exposure area to fresh air.

**After skin contact:** Immediately wash skin with soap and water for at least 20 minutes while removing contaminated clothing.

**After eye contact:** Immediately flush with cool water for at least 20 minutes.

**After ingestion:** Do NOT induce vomiting. Dilute by rinsing the mouth and giving cool water to drink. Never give anything by mouth to an unconscious person

### SECTION 5: FIREFIGHTING MEASURES

#### Suitable (and unsuitable) extinguishing media

**Suitable extinguishing media:** Use extinguishing agent suitable for surrounding fire.

**Unsuitable extinguishing media:** None.

**Specific hazards arising from the chemical:** Not flammable or combustible.

#### Special protective equipment and precautions for firefighters

In case of fire, wear a full-face positive-pressure self-contained breathing apparatus and protective suit. Shut off all ignition sources and apply cooling water to sides of containers that are exposed to flames until well after fire is out. Avoid breathing vapors, keep upwind. Keep unnecessary personnel away.

**Neutral Non-enzymatic Pretreatment Spray****SECTION 6: ACCIDENTAL RELEASE MEASURES****Personal precautions, protective equipment, and emergency procedures**

Wear skin, eye, and respiratory protection. Stop or reduce any leaks if it is safe to do so. Spills present a slipping hazard. Keep unnecessary personnel away. Ensure clean-up is conducted by trained personnel only. Make sure area is slip-free before re-opening to traffic.

**Environmental Precautions**

Avoid runoff into sewers, ditches, and waterways.

**Methods and materials for containment and cleaning up**

Contain and absorb with sand or other absorbent material and place into containers for later disposal. Wash site of spillage thoroughly with water. In case of large spills, dike far ahead of spill to prevent further movement. Recover by pumping or by using a suitable absorbent material and place into containers for later disposal. Dispose in suitable waste container.

**SECTION 7: HANDLING AND STORAGE****Precautions for safe handling**

Use with adequate ventilation. Avoid breathing vapors, spray, or mist.

**Conditions for safe storage, including any incompatibilities**

Store away from incompatible substances and avoid extreme temperatures. Keep the containers tightly closed when not in use.

**SECTION 8: EXPOSURE CONTROLS / PERSONAL PROTECTION**

OSHA permissible exposure limit (PEL), American Conference of Governmental Industrial Hygienists (ACGIH) Threshold Limit Value (TLV), and any other exposure limit used or recommended by the chemical manufacturer, importer, or employer preparing the safety data sheet, where available.

**Appropriate engineering controls**

Good general ventilation should be sufficient. Handle in accordance with good industrial hygiene and safety practices.

**Individual protection measures**

**Respiratory:** If respiratory protection is required, it must be based on the contamination levels found in the workplace and be jointly approved by the National Institute for Occupational Safety and Health and the Mine Safety and Health Administration (NIOSH-MSHA) and must not exceed the working limits of the respirator.

**Eye:** Safety glasses are recommended but not required.

**Hand/Body:** Gloves are recommended but not required.

**SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES**

Appearance	: Clear, light straw
Odor	: Mild
Odor threshold	: Not available
pH	: 6.5 - 8.5 (neat)
pH	: Not available
Melting point / range	: Not available
Boiling point / range	: Not determined
Flash point	: Not determined
Evaporation rate	: Not determined
Flammable/Explosive limits - lower	: Not available
Flammable/Explosive limits - upper	: Not available

## Safety Data Sheet

### Neutral Non-enzymatic Pretreatment Spray

Vapor pressure	: Not determined
Vapor density	: Not determined
Solubility in water	: Completely soluble
Partition coefficient (n-octanol/water)	: Not available
Autoignition temperature	: Not available
Decomposition temperature	: Not available
Viscosity	: Not available
VOC content	: Not available
Density	: Not available
Relative Density	: Approximately 1.015 (Specific Gravity)

## SECTION 10: STABILITY AND REACTIVITY

<b>Reactivity</b>	: None known.
<b>Chemical stability</b>	: Stable under normal conditions.
<b>Possibility of hazardous reactions</b>	: None known.
<b>Conditions to avoid</b>	: Avoid extremes of temperature.
<b>Incompatible materials</b>	: None known.
<b>Hazardous decomposition products</b>	: Not available.

## SECTION 11: TOXICOLOGICAL INFORMATION

### Likely routes of exposure including symptoms related to characteristics

**Inhalation:** May cause allergic respiratory reaction.

**Skin contact:** Health injuries are not known or expected.

**Eye contact:** Health injuries are not known or expected.

**Ingestion:** Ingestion of large quantities may cause gastrointestinal irritation with nausea, vomiting, and diarrhea.

**Physical/Chemical:** No physical/chemical hazards are anticipated for this product.

### Other relevant toxicity information:

The use of this product by the end-user is safe under normal and reasonable foreseen use.

### Numerical measures of toxicity, including delayed and immediate effect

#### Carcinogenicity information

Hazardous Component(s)	NTP Carcinogen	IARC Carcinogen	OSHA Carcinogen
Not known	Not known	Not known	Not known

**Carcinogenicity** : None of the ingredients in this product are listed as carcinogens by the International Agency for Research on Cancer (IARC), the National Toxicology Program (NTP), or the Occupational Safety and Health Administration (OSHA).

**Mutagenicity** : None of the ingredients in this product are known to cause mutagenicity.

**Reproductive Toxicity** : Not available.

# Safety Data Sheet

## Neutral Non-enzymatic Pretreatment Spray

### SECTION 12: ECOLOGICAL INFORMATION

#### Aquatic Toxicity:

The following toxicity information is available for the hazardous ingredient(s) when used as technical grade. The product is anticipated to be safe for the environment at concentrations predicted in settings under normal-use conditions. One component, sodium tetraborate pentahydrate, is considered harmful to the aquatic environment when used in high concentrations.

#### Toxicity to fish:

The product is not anticipated to be toxic to fish. The following data is for the product tested in accordance with EPA-821-R-02-012 (Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms).

Hazardous substance	Parameter	Value	Exposure Time	Species
Product	LC50	16.21 mg/L	96 hrs.	Pimephales promelas

#### Toxicity to aquatic invertebrates:

The product is not anticipated to be toxic to aquatic invertebrates. The following data is for the product tested in accordance with EPA-821-R-02-012 (Methods for Measuring the Acute Toxicity of Effluents and Receiving Waters to Freshwater and Marine Organisms).

Hazardous substance	Parameter	Value	Exposure Time	Species
Product	LC50	57.43 mg/L	96 hrs.	Ceriodaphnia dubia

#### Toxicity to algae:

The product is not anticipated to be toxic to algae.

#### Persistence and degradability:

The majority of the components of this product are anticipated to be biodegradable.

#### Bioaccumulation Potential:

The bioaccumulation potential of this product has not been determined.

#### Mobility:

The mobility of this product (in soil and in water) has not been determined.

### SECTION 13: DISPOSAL CONSIDERATIONS

- Waste Number and Description** : Not applicable
- Disposal Considerations** : Not regulated
- Disposal of products** : This product is not a RCRA hazardous waste and can be disposed of in accordance with federal, state, and local regulations.
- Disposal of packages** : Do not reuse container. Rinse it with water and place it in trash.

### SECTION 14: TRANSPORT INFORMATION

The information in this section is for reference only and should not take the place of a shipping paper (bill of lading) specific to an order. Please note that the proper shipping classification may vary by packaging, properties, and mode of transportation.

- Land Transport (DOT) : Not dangerous goods
- Sea Transport (IMDG/IMO) : Not dangerous goods
- Air transport (IATA) : Not dangerous goods

# Safety Data Sheet

## Neutral Non-enzymatic Pretreatment Spray

### SECTION 15: REGULATORY INFORMATION

#### Occupational Safety and Health Act:

Hazard Communication Rule, 29 CFR 1910.1200: The Occupational Safety and Health Administration (OSHA) require Safety Data Sheets (SDSs) to provide information about any hazard that may be associated with the product and make this information available in the workplace. Since the use pattern and exposure in the workplace are generally not consistent with those experienced by consumers, this SDS may contain health hazard information not relevant to consumer use.

#### United States Regulatory Information

**TSCA 12 (b) Export Notification** : All components are listed.  
**CERCLA/SARA Section 302 EHS** : None  
**CERCLA/SARA Section 311/312** : Acute health hazard  
**CERCLA/SARA Section 313** : Does not contain any components with reportable quantities.  
**California Proposition 65** : This product does not contain substances listed under California Proposition 65.

**California Cleaning Product Act to Know Act of 2017 (SB 258)** : This product is labeled in compliance with CA SB258.

Chemical Name	CAS Registry Number	Functional Purpose	List(s)
Water	7732-18-5	Diluent	Not Applicable
Benzotriazole	95-14-7	Corrosion Inhibitor	Not Applicable
Octylphenoxy poly(ethyleneoxy)ethanol	9002-93-1	Cleaning Agent	Not Applicable
Glycerin	56-81-5	Stabilizer	Not Applicable
Sodium Xylene Sulfonate	1300-72-7	Cleaning Agent	Not Applicable
Tetrasodium EDTA	64-02-8	Chelating Agent	Not Applicable
Cocamidopropyl Betaine	61789-40-0	Cleaning Agent	Not Applicable
DMDM Hydantoin	6440-58-0	Biocide	Not Applicable
Colorant	Not Available	Dye	Not Applicable

\*refer to [www.belimed.com/en-us/solutions/consumables-protect](http://www.belimed.com/en-us/solutions/consumables-protect) for electronic links to designated lists

#### The ingredients of this product are reported on the following inventories:

#### TSCA 8 (b) Inventory Status

All components of this product are either listed on or exempt from the U.S. Toxic Substances Control Act (TSCA) chemical substance inventory.

### SECTION 16: ADDITIONAL INFORMATION

**NFPA health hazard:** 1

**NFPA fire hazard:** 0

**NFPA reactivity:** 0

**Additional information:** None

**DISCLAIMER:** The information provided on this sheet is not to be considered as warranty or quality specification. The information is given to provide general knowledge as to health and safety based on our understanding of safe handling, storage, transportation, disposal, and release. The information relates to the specific material designated and may not be valid for such material used in combination with any other materials or any process unless specified in the text.

**This safety data sheet contains changes from the previous version in sections:** company logo

**Prepared by:** Belimed Inc.

**Issue date:** August 7, 2024