

Product Specification

Closed Case Carts

Product Description

Belimed closed case carts are specifically designed to be an integral part of the instrument transport process. The line includes 6 cart sizes, with a wide range of available heights, widths, and lengths.

Application

For use in healthcare facilities to facilitate secure transport of sterile instrument sets from the SPD to the OR, and hygienic return of contaminated sets back to the SPD for cleaning, decontamination, and sterilization.

Standard Features

Construction/Design

These case carts employ seam-welded cabinet construction, with an internal structural bracing skeleton that provides additional rigidity and support.

- Type 304 stainless steel construction, 16 gauge, #4 finish
- Interlocked internal frame construction
- Non-marring wrap-around bumper
- Positive-lock door latches
- Heavy-duty door tie-back latches
- Integrated push handles for easy maneuvering
- Fully insulated doors and top
- 1 or 2 standard 18 gauge perforated adjustable stainless shelves (quantity depends on model)
- 8.5" x 11" card/documentation holder
- 6" diameter all-stainless steel casters, with sealed bearings
- Tamper evidence security tab

Dimensions

Units are measured by case (internal) and outer (total) dimensions:

Model	Case Dimensions (inches)			Outer Dimensions (inches)			Doors
	Width	Depth	Height	Width	Depth	Height	
MMCCB3	40.25	25.00	32.00	45.75	26.88	39.63	2
MMCCB4	40.25	25.00	47.00	45.75	26.88	54.63	2
MMCCB6	23.63	25.00	32.00	29.13	26.88	39.63	1
MMCCB9	23.63	25.00	47.00	29.13	26.88	54.63	1
MMCCB12	36.00	25.00	32.00	41.50	26.88	39.63	2
MMCCB14	46.50	25.88	32.00	52.00	27.75	39.63	2

Outer dimensions include integral wrap-around bumper base.



Benefits

Stability & Durability

Rigid heavy-gauge stainless steel construction with seam-welded components and full interlocked skeleton frame assures years of reliable operation without distortion, cracks, or other damage. Cart-washer compliant stainless steel sealed-bearing casters are standard.

Capacity that fits your facility

Since not all set containers and container quantity requirements are the same, we offer carts with the widest range of internal dimensions. Minimum depth carts for space-constrained facilities, and maximum depth carts to accommodate cases with high set count. Wrap-friendly solid stainless steel adjustable shelves are standard on all units.

Features you need are standard

Positive-lock door latches keep contents secure, and push/pull door tiebacks make loading/unloading easy, and assure complete unit washing and drying. Documentation card holder, shelves, and uprated casters are all standard, not costly options.



For more information, visit belimed.com

Closed Case Carts

Accessories

Glass Doors

These allow you to check the contents of your cart with just a glance. The high impact polycarbonate glass features high temperature resistance along with high clarity.

Shelves

- Solid Adjustable Roll Out Shelf
- Additional adjustable perforated
- Roll-out full-extension adjustable perforated
- Flip-up top extensions

Vertical Handles

- Cool Touch option for Vertical handles

Casters

- Stainless steel with directional and rolling locks

Security Aids

- Perimeter silicone door gaskets
- Door security tabs
- Flag system
- Clean Dirty Indicator
- Tow Hitch



Distributed by Belimed for Mac Medical

REV7.6.2023_Q0

Engineers of Confidence.

Belimed
Infection Control

Belimed, Inc.
1.800.451.4118 | 1.843.216.7424
belimed.com