



WD 290 Washer-Disinfector

The WD 290 optimizes resources and increases workflow efficiency. Making it the ideal choice in planning for tomorrow's needs and setting higher expectations for your CSSD.

Engineers of Confidence.



WD 290 Washer–Disinfecter

Safety Features

The WD 290 reaches a thermal disinfection level of A_0 3000 Standard versus the typical disinfection level of A_0 600. In order to minimize the possibility of cross contamination, sensors in the system will not release the load unless all cycle parameters are met.

Workflow Efficient

The high performance and large capacity of the washer ensure that instruments are getting to the clean side at just the right time, balancing workflow through the CSSD. The longer chamber allows processing of larger orthopedic and bariatric sets. Together, these features eliminate bottlenecks at prep and pack and streamline the department.

Resource Saving

With a combination of Belimed's SmartFill technology, which senses the load size and adjusts water use accordingly, and a soaking wash approach, the WD 290 offers the lowest water use per tray, in the industry. Innovative design makes the WD 290 completely serviceable from the front and back of the unit; no side service area required. This results in a smaller width and overall footprint, conserving another valuable resource in the CSSD: Space.

Configurations & Options

- Standard Configuration:
 - Drain water cool down
 - Three (3) detergent pumps and flow meters
 - Process status display for clear visibility throughout the department
- Steam or Electric Heating
- Manual Loading



Accessories

Belimed offers standard instrument racks of various levels, along with specialty racks, to accommodate all instrument reprocessing needs. The most commonly used rack is the 3-level instrument rack. Each level can hold up to 3 instrument trays. Transport carts are also available to easily move racks around the department.



Instrument Rack	Tray Capacity	Useable Space
1-Level	3	Bottom level: 23.6"
2-Level	6	Each level: 10.6"
3-Level	9	Each level: 6.3"
4-Level	12	Each level: 4.1"
5-Level	15	Each level: 2.7"

REV2.7.2022_Q320

Engineers of Confidence.