



Belimed
Infection Control

CS 750 Cart Washer

Resource responsible, space saving, and workflow efficient. These define the Belimed CS 750 cart washer. It's the only real choice to meet today's demand for fast, effective and efficient cart and container processing in the CSSD.

Engineers of Confidence.



CS 750 Cart Washer

Workflow Efficient

Our high output dual steam heated dryers, tilting floor and high temperature rinse ensure that drying time is reduced to an absolute minimum. Carts are cleaned and dried in approximately 8 minutes so the clean side staff can focus on the instruments instead of hand drying wet carts. The large, bright status indicator, easily visible from across the room, lets personnel know that the cycle is complete and ready to be unloaded, further improving productivity. Let your washer/disinfectors take care of instruments while the CS 750, optimized for washing and drying as quickly as possible, handles your carts.

Space Saving

Tinted glass pivot doors, compact wash tank and plumbing design increases chamber volume while reducing the overall footprint. As a result, valuable floor space is left available for additional needs of the department.

Resource Responsible

A highly efficient tank wash system, coupled with an on-demand rinse water heater, reduces water use to about 8 gallons per cycle. This unique combination also eliminates the need to cool drain discharge water, conserving water and reducing operating costs.

Accessories

Carts for sterilization containers and lids allow you to get more out of your cart washer and free up the washer/disinfectors for instrument loads. Improved workflow, more efficient.

Container Cart Series

- Stainless steel casters with non-marking 5" diameter neoprene wheels
- Dimensions of cart: 26"W x 49"L x 67"H
- 5 container levels
- Bottom level for lids, with 6 or 12 dividers
- 20 container capacity (7" tall containers)
- Shelves are easily removable for taller containers



Tilting floor enables water runoff and faster cart drying.

Parameter	CS 750 S	CS 750 M	CS 750 L
Usable chamber dimensions H x W x D, in (mm)	78.7" x 39.4" x 53.1" (2000 x 1000 x 1350)	78.7" x 39.4" x 90.6" (2000 x 1000 x 2300)	78.7" x 39.4" x 116.1" (2000 x 1000 x 2950)
Outer dimensions H x W x D, in (mm)	114.2" x 120" x 88.6" (2900 x 3050 x 2250)	114.2" x 120" x 126.0" (2900 x 3050 x 3200)	114.2" x 120" x 151.6" (2900 x 3050 x 3850)
Pit depth, in (mm)	4.7" (120)	4.7" (120)	4.7" (120)

REV2.4.2022_Q0

Engineers of Confidence.