

Product Description

Belimed Protect® Pretreatment Foamer

Product Description

The Belimed Protect® Pretreatment Foamer is designed to easily and ergonomically apply Belimed Protect Pretreatment to surgical instruments. The foamer includes a unique cart for mobility and two wands, one for general trays and one for cannulated instruments.

Standard Features


- The Pretreatment Foamer allows for ergonomic application for staff, resulting in better compliance for pretreating instruments
- Two Application Wands – one for general instruments and one specially designed for cannulated instruments
- Triple the coverage of manual spraying in the same application time
- The foam created from the Pretreatment Foamer will not aerosolize and is safe to use in the Operating Room
- The cart design holds 1 gallon or 2.5 gallon container sizes, which is more cost effective and environmentally friendly than traditional manual spray bottles

Specifications

- Part Number: 1064361
- Working Voltage: 110-120VAC
- Power Consumption: 30W
- Frequency: 47-63Hz



Item Number	Description
QFC-I-01	Power Supply
QFC-I-02	Fuse
QFC-I-03	Pressure Control
QFC-I-04	Air Pump
QFC-I-05	Foaming Out Port
QFC-I-06	Silencer
QFC-I-07	Control Box
QFC-I-08	Foam Generator
QFC-I-09	Water Pump
QFC-I-10	Detergent Entry
QFC-I-11	Check Valve
QFC-O-01	Cabinet
QFC-O-02	Power Outlet
QFC-O-03	Hose Connection
QFC-O-04	Power Switch
QFC-O-05	Liquid Wash



For more information, visit belimed.com

DEV12.15.2021_Q0

Engineers of Confidence.

Belimed
Infection Control

Belimed, Inc.
1.800.451.4118 | 1.843.216.7424
belimed.com