



Belimed Protect® Pretreatment Foamer

The Belimed Protect® Pretreatment Foamer is designed to easily and ergonomically apply Belimed Protect Pretreatment to surgical instruments. The foamer includes a unique cart for mobility and two wands, one for general trays and one for cannulated instruments.

Engineers of Confidence.



Belimed Protect® Pretreatment Foamer

Standard Features

- The Pretreatment Foamer allows for ergonomic application for staff, resulting in better compliance for pretreating instruments
- Two Application Wands – one for general instruments and one specially designed for cannulated instruments
- Triple the coverage of manual spraying in the same application time
- The foam created from the Pretreatment Foamer will not aerosolize and is safe to use in the Operating Room
- The cart design holds 1 gallon or 2.5 gallon container sizes, which is more cost effective and environmentally friendly than traditional manual spray bottles

Specifications

- Part Number: 1064361
- Working Voltage: 110-120VAC
- Power Consumption: 30W
- Frequency: 47-63Hz



Easy Access Storage



Solution for Cannulated Instruments



Ergonomic Application



Triple the Coverage

DEV12.15.2021_Q0

Engineers of Confidence.