



Warming Cabinets

User friendly, durable, and versatile. Belimed's warming cabinets are precisely engineered to be simple to operate and maintain while providing years of worry-free material storage. Our multiple options and sizes fit any cabinet need for safe, consistent warming of medical liquids and linens.

Engineers of Confidence.



Warming Cabinets

User Friendly

Easy to use touch panel controls allow the user to easily modify temperature set points but a key lockout feature prevents unauthorized alterations. The audible/visual over-temperature alarm minimizes the possibility of overheated supplies being delivered to the surgical suite. The "Easy Log" datalogger lets users to conveniently download temperature data recorded as far back as 2 years.

Durable

Sturdy cabinet construction and durable components ensure years of service life. Easy access electrical components greatly simplify system maintenance and periodic upkeep.

Versatile

Multi-chamber models feature independent access doors and individual temperature controls for each chamber to allow storage of dissimilar materials in the same warmer.

Distributed by Belimed for Mac Medical

Features

A variety of options and accessories are available to enhance our warmers' capabilities and workflow value:



Standard

- Stainless steel cabinet and shelves
- Choice of glass or stainless steel doors
- Audible and visual over-temperature alarm, and keyed door/temperature lockout
- Integrated data logger with removable USB
- Fully Insulated for uniform heating
- Independent digitally controlled chambers

Options

- Pass-through, Recessed, or Standalone configurations
- Left hinge with right door swing
- Sloping top
- Electronic keypad lock
- Antimicrobial copper alloy handle
- Premium version with Intuitive touch software and Ethernet/WiFi capability for remote access

Accessories

- Additional Stainless steel perforated shelves
- Mobile base - 2 and 4 inches
- Mobile stand (Tabletop models only)
- Mobile cabinet stand (Tabletop models only)
- Roll-out basket
- Seismic anchoring kits



Specifications

Model	Chamber Type	Storage Capacity (cu. ft.)	Exterior Dimensions, W x D x H (inches)	Shelves
SWC1518	Single	1.3	18.0 x 17.25 x 22.0	None std, 1 optional
SWC1824	Single	3.0	24.0 x 20.5 x 24.5	None std, 1 optional
SWC24	Single	5.27	30.0 x 26.5 x 24.5	1, adjustable
SWC36	Single	8.65	30.0 x 26.5 x 36.0	2, adjustable
SWC2460	Single	9.74	24.0 x 20.5 x 64.5	3, adjustable
SWC60	Single	13.05	30.0 x 20.5 x 64.5	3, adjustable
SWC72	Single	21.1	30.0 x 26.5 x 74.5	2, adjustable, 1, Fixed
DWC2460	Double	7.47 (2.65 upper/4.82 lower)	24.0 x 20.5 x 64.5	3, adjustable
DWC60	Double	9.72 (3.45 upper/6.27 lower)	30.0 x 20.5 x 64.5	3, adjustable
DWC2424	Double	13.24 (4.06 upper/9.18 lower)	24.0 x 26.5 x 74.5	3, adjustable
DWC24	Double	17.21 (5.28 upper/11.93 lower)	30.0 x 26.5 x 74.5	3, adjustable
DWC36	Double	17.64 (8.65 upper/8.99 lower)	30.0 x 26.5 x 74.5	4, adjustable
TWC24	Triple	14.01 (3.89 upper/3.72 middle/6.4 lower)	30.0 x 26.5 x 78.2	1, adjustable

REV12.21.2020_Q320

Engineers of Confidence.