



Cleaning Indicator

The Belimed Protect Cleaning Indicator (BPCI) is a verification test used to monitor the effectiveness of the cleaning process in washer / disinfectors. The BPCI and special holder mimic soiled surgical instruments, providing an appropriate challenge for all factors in automated cleaning including time, temperature, pressure, detergent, and water quality.

Engineers of Confidence.



Belimed Protect® Cleaning Indicator

The Belimed Protect Cleaning Indicator (BPCI) is a verification test used to monitor the effectiveness of the cleaning process in washer / disinfectors. The BPCI and special holder mimic soiled surgical instruments, providing an appropriate challenge for all factors in automated cleaning including time, temperature, pressure, detergent, and water quality.

Product Features


- The Belimed Protect Cleaning Indicator was thoroughly tested to ensure the coupon and holder provide a sufficient challenge and consistent monitoring results.
- Can be clipped to the side of surgical trays during testing, providing a better challenge for blind spots in the washer disinfectant.
- BPCI is made of a synthetic recyclable substrate which is more sustainable and cost-effective than metal-based indicators.
- BPCI has a three-year shelf life, more than double the shelf life of other cleaning verification products.

CLEANING INDICATOR FOR CLEANING PERFORMANCE TESTS

Designed for routine control of cleaning processes and for monitoring automatic washers

HOLDER

Designed to mimic the cleaning challenges posed by a large process load during the washing cycle



Part Number	Description
1072811	BPCI Cleaning Indicator Card - Pouch of 200
1072812	BPCIH Protect Cleaning Indicator Holder



Developed by  Terragene®

DEV10.14.2022_Q0

Engineers of Confidence.