



## Container Cart

The CWA-004 Container Cart provides an effective method for cleaning containers with a cart washer. Utilizing a cart washer for container washing increases overall container throughput and reduces the demand on instrument washer / disinfectors.

Engineers of Confidence.



# Container Cart

The CWA-004 Container Cart provides an effective method for cleaning containers with a cart washer. Utilizing a cart washer for container washing increases overall container throughput and reduces the demand on instrument washer / disinfectors.

**Features:**

- Manufactured from stainless steel for long life
- Donut bumper on each leg
- Stainless steel casters with non-marking 5" diameter neoprene wheels
- Fully welded frame
- Dimensions of cart: 26" W x 49" L x 67" H
- 4 levels
- Capacity for 16 containers (picture does not include 4th container shelf)
- Includes an additional 6 lid dividers for a 12 lid capacity

Part Number	Dimensions
CWA-004	26" W x 49" L x 67" H

Distributed by Belimed for Pedrigo



REV10.28.2022\_Q0

Engineers of Confidence.