

Preventive Maintenance Kits (PM Kits)

To order PM Kits for regular delivery, please contact Belimed's Parts Admin Team at 1 (800) 451-4118, Option 2 or via e-mail at parts.us@belimed.com

Washer / Disinfectors

Part Number	PM Kit Description	Maintenance Cycle Time
107-4145	PM Pack, WD200, Annual	Every 3,000 Cycles
107-4146	PM Pack, WD200, Quarterly	Every 1,500 Cycles
107-4141	PM Pack, WD230 250, Plastic Hub, Annual	Every 3,000 Cycles
107-4142	PM Pack, WD230 250, Plastic Hub, Quarterly	Every 1,500 Cycles
107-4143	PM Pack, WD230 250, SS Hub, Annual	Every 3,000 Cycles
107-4144	PM Pack, WD230 250, SS Hub, Quarterly	Every 1,500 Cycles
106-9944	PM Pack, WD250-12	Every 3,000 Cycles
107-4138	PM Pack, WD290, Annual	Every 3,000 Cycles
107-4140	PM Pack, WD290, Quarterly	Every 1,500 Cycles
1052056	PM Pack, WD290IQ Annual	Every 3,000 Cycles
1052057	PM Pack, WD290IQ Quarterly	Every 1,500 Cycles

Cart Washers

Part Number	PM Kit Description	Maintenance Cycle Time
107-4147	PM Pack, CS720 ES387, Annual with Seals	*Not Cycle Based
107-4148	PM Pack, CS720 ES387, Annual Incl. Hoses	*Not Cycle Based
107-4149	PM Pack, CS720 ES387, Quarterly	*Not Cycle Based
107-4150	PM Pack, CS721, Annual with Seals	*Not Cycle Based
107-4151	PM Pack, CS721, Annual Incl. Hoses	*Not Cycle Based
107-4152	PM Pack, CS721, Quarterly	*Not Cycle Based
107-4153	PM Pack, CS750, Annual with Seals	Every 6,000 Cycles
107-4154	PM Pack, CS750, Annual Incl. Hoses	Every 3,000 Cycles
107-4155	PM Pack, CS750, Quarterly	Every 1,000 Cycles
1039696	PM Pack, CS750WAY, Quarterly	Every 1,500 Cycles
1039697	PM Pack, CS750WAY, Annual	Every 3,000 Cycles
1025411	PM Pack, WD750 Quarterly	Every 1,000 Cycles
1025432	PM Pack, WD750, Annual	Every 3,000 Cycles

Sterilizers

Part Number	PM Kit Description	Maintenance Cycle Time
107-4865	PM Pack, 446 VS1, Annual	Every 2,600 Cycles
107-4866	PM Pack, 446 VS1, Quarterly	Every 650 Cycles
107-4867	PM Pack, 446 VS1, Semi-Annual	Every 1,300 Cycles
107-4851	PM Pack, 559VS1, Annual	Every 2,600 Cycles
107-4852	PM Pack, 559VSx, Quarterly	Every 650 Cycles
107-4853	PM Pack, 559VS2, Annual	Every 2,600 Cycles
107-4854	PM Pack, 66x VS1, Annual	Every 2,600 Cycles
107-4855	PM Pack, 66x VSx, Quarterly	Every 650 Cycles
107-4856	PM Pack, 66x VS2, Annual	Every 2,600 Cycles
107-4857	PM Pack, 96xx HS1, Annual	Every 2,600 Cycles
107-4858	PM Pack, 96xx HSx, Quarterly	Every 650 Cycles
107-4859	PM Pack, 96xx HS2, Annual	Every 2,600 Cycles
107-4860	PM Pack, GR96xx HS1, Annual	Every 2,600 Cycles
107-4861	PM Pack, GR96xx HSx, Quarterly	Every 650 Cycles
107-4862	PM Pack, GR96xx HS2, Annual	Every 2,600 Cycles
107-4863	PM Pack, MST-V 336, Annual	Every 2,600 Cycles
107-4864	PM Pack, MST-V 336, Quarterly	Every 650 Cycles
109-8961	PM Pack, MSTH10 96xx HS2, Annual	Every 2,600 Cycles
109-8962	PM Pack, MSTH10 96xx HS1, Annual	Every 2,600 Cycles
109-8963	PM Pack, MSTH10 96xx HSx, Quarterly	Every 650 Cycles

Accessories

Part Number	PM Kit Description	Maintenance Cycle Time
1019560	PM Pack, Knight Sink Pump, Annual	*Not Cycle Based
1045521	PM Pack, ELD 16 30 and 45kW, Annual	Every 2,600 Cycles

*Refer to the Equipment Manual for specific maintenance schedule.

REV4.11.2023_Q0

Engineers of Confidence.