

PREP & PACK WORKSTATION - TECHNICAL SPECIFICATIONS

REQUIRED INPUTS

Base Type	Description	Value	Utility Notes	Connection Method	Final Utility Conn. B	
					Others	Belimed
CA (OPTIONAL)	Utility:	Compressed Air	12' Polyurethane retractable (coiled) hose with 1/4" ID provided by Belimed. Air must be dry and oil free. If configured with air tool kit.	Flex		X
	Connection Size/Type:	1/4 inch / MNPT				
	Pressure (Max):	30 psig				
EL (OPTIONAL)	Utility:	Electrical - Base	Powered base frame, includes 9.5 ft. line cord with NEMA 5-15 plug, center exit 2" above finished floor.	Receptacle		X
	Voltage:	120 V				
	Cycles:	60 Hz				
	Current/Watts:	1.6 A / 192W				

OPTIONAL OUTLET POWER BAR FOR ALL BASE TYPES

	Utility:	Electrical	8 outlet power bar, includes 8 ft. line cord with NEMA 5-15 plug, right exit.
	Voltage:	120 V	
	Cycles:	60 Hz	
	Current/Watts:	15 A / 1800W	

PHYSICAL SPECIFICATIONS

FIXED	30" X 48":	160 lbs*	Maximum static load capacity 1000 lbs.
	30" X 60":	189 lbs*	
	30" X 72":	220 lbs*	
	36" X 48":	177 lbs*	
	36" X 60":	215 lbs*	
	36" X 72":	257 lbs*	

PHYSICAL SPECIFICATIONS

ADJUSTABLE	30" X 48":	192 lbs*	Maximum static load capacity 1000 lbs.
	30" X 60":	228 lbs*	
	30" X 72":	261 lbs*	
	36" X 48":	207 lbs*	
	36" X 60":	252 lbs*	
	36" X 72":	287 lbs*	

PHYSICAL SPECIFICATIONS

POWERED	30" X 48":	218 lbs*	Maximum frame lift capacity 500 lbs.
	30" X 60":	254 lbs*	
	30" X 72":	291 lbs*	
	36" X 48":	223 lbs*	
	36" X 60":	272 lbs*	
	36" X 72":	314 lbs*	

*Weights calculated with no accessories.

Additional Information: